

SPO2 Analyzer-Testers

OXSIM



Compact Portable SPO2 Simulator

The OxSim is a optical, pocket sized unit that generates four saturations, three pulse rates and a low perfusion mode all at the touch of a button. This makes it an indispensable tool for any Biomed.

[See More...](#)

SMARTSTAT



Comprehensive SPO2 Analyzer

SmartSat is the only universal SpO2 simulator which uses the manufacturer-approved "direct electrical simulation" method to provide unmatched accuracy and reliability. In fact, SmartSat is used by leading pulse oximeter manufacturers and is emerging as the industry standard pulse oximetry tester.

FINGERSIM



Dye Induced Reference Phantoms

Start Up Kit-3 FingerSims Includes the Following:

80% FingerSim

90% FingerSim

97% FingerSim

FingerSim Holder

Plastic Carrying Case (for FingerSims, User Manual & FingerSim Holder)

User Manual



OxSim - An Optical Pocket SPO2 simulator that runs on one AA alkaline battery or drops into the **SimCube** carrying case and ties directly into the **SimCube** Battery Boost Option

SPO2 Simulations will never be the same!!!!!!

- ◆ Optical oxygen saturation simulations
- ◆ Automatic detection of sensor type
- ◆ 4 different Saturation rates
- ◆ 3 different Pulse Rates
- ◆ Low Perfusion Simulation
- ◆ IR and Red LED detection
- ◆ Runs on one AA alkaline battery
- ◆ Drop tested 40 times with no failures